

SUBSTANCE	SCREENING THRESHOLD	CONFIRMATION THRESHOLD	DETECTION TIME IN URINE (estimate only)
AMPHETAMINE (AMP)	500 ng/mL (2000 ng/mL= High)	100 ng/mL	2 to 3 days
BARBITURATES (BAR)	200 ng/mL (1000 ng/mL=High)	50 ng/mL	1 day to 3 weeks
BENZODIAZEPINES (BZO)	200 ng/mL (1000 ng/mL= High)	25 ng/mL	Up to 7 days
BUPRENORPHINE II (BUP)	10 ng/mL (>100 ng/mL=High)	5 ng/mL Buprenorphine 10 ng/mL Norbuprenorphine	2 to 4 days
CANNABINOID (THC)	50 ng/mL (200 ng/mL=High)	5 ng/mL	Up to 6 weeks (regular/chronic use) 2 to 3 days (occasional use)
COCAINE (COC)	150 ng/mL (1000 ng/mL=High)	100 ng/mL	2 to 4 days
ETHYL GLUCURONIDE (EtG)	500 ng/mL (2000 ng/mL=High)	250 ng/mL Ethyl Glucuronide 100 ng/mL Ethyl Sulfate	Up to 80 hours
ETHANOL (ETOH)	20 mg/dL (600 mg/dL=High)	20 mg/dL	Dose dependent
FENTANYL (FEN)	2 ng/mL	2.5 ng/mL Fentanyl, Acetylfentanyl 5 ng/mL Norfentanyl	1 to 3 days
HEROIN (6-ACETYLMORPHINE)	10 ng/mL	10 ng/mL	24 hours
HYDROCODONE	300 ng/mL	25 ng/mL Hydrocodone, Hydromorphone 50 ng/mL Norhydrocodone	2 to 4 days
ECSTASY (MDMA)	500 ng/mL (1000 ng/mL=High)	100 ng/mL	1 to 2 days
METHADONE (MTD)	300 ng/mL (1000 ng/mL=High)	25 ng/mL	7 to 10 days

OPIATES (OPI)	300 ng/mL (1000 ng/mL-High)	25 ng/mL Codeine, Morphine	2 to 4 days
OXYCODONE (OXY)	100 ng/mL 1000 ng/mL=High)	25 ng/mL Oxycodone, Oxymorphone; 50 ng/mL Noroxycodone	2 to 4 days
PHENCYCLIDINE (PCP)	25 ng/mL (100 ng/mL-High)	25 ng/mL	Up to 10 days
SYNTHETIC MARIJUANA K2/SPICE	10 ng/mL	1 ng/mL	Up to 72 hours
GABAPENTIN (GABA)	1,000 ng/mL	500 ng/mL	1 to 2 days
DEXTROMETHORPHAN (DXM)	Screen not available	2.5 ng/mL	2 to 4 days
TRAMADOL (TRA)	100 ng/mL	100 ng/mL	2 to 4 days
BATH SALTS (BATH)	Screen not available	25 ng/mL	2 to 3 days
KRATOM (KRA)	Screen not available	1 ng/mL	Up to 3 days
KETAMINE (KET)	Screen not available	25 ng/mL	1-4 days

ZOLPIDEM	20 ng/mL	20 ng/mL	1-3 days
----------	----------	----------	----------